## Applause® – Directions For Use

GB

### Indications: Dental Ceramic Alloy, Type 4

**Contraindications:** Do not use in patients with known hypersensitivity to one or several metals contained in the alloy.

**Warnings/Precautions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Use proper safety equipment and a certified industrially ventilated work area when processing this alloy.

**Adverse Reactions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

### Instructions:

**Investment:** Phosphate.

**Burn out:** 700 °C (1300 °F), 30–120 min.

**Cast at:** 1345 °C (2453 °F), do not use carbon flux or crucibles.

**Finish with:** Carbide burs or ceramic bound stones.

**Pre-Clean:** Blast with alumina, 50–125 µm/< 3 bar (44 psi), steam clean or ultrasonic with distilled water.

**Oxidation:** 650 °C (1200 °F) to 1038 °C (1900 °F); full vacuum, 10 min. hold.

**Post-Clean:** Blast with alumina, 100–125 µm/< 2 bar (30 psi), steam clean or ultrasonic with distilled water.

**Heat treatment:** None.

**Recommended Porcelain:** Ceramco® II Silver, Duceram® Kiss, Duceram® Plus.

**Pre-Solder:** SMG®-YW/Ceramco® White; Flux: Floride.

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