MSDS

Rx only

**DENTSPLY**  
CERAMCO

# Option®

Dental Ceramic Alloy, Type 4/Alliage céramique dentaire, Type 4/  
Leghe ceramiche dentale, Tipo 4/Dental-Keramiklegierung, Typ 4/  
Aleación cerámica dental, Tipo 4

2.0 % Au, 78.8 % Pd

CTE/CET/WAK 14.3 (500 °C), 14.6 (600 °C)

**REF 1813 0005****1 T.oz/31.1 g**

DENTSPLY CERAMCO  
DENTSPLY International Inc.  
570 West College Ave.  
York, PA 17405-0872  
1-800-487-0100  
1-717-699-4190

Color/Couleur/Colore/Farbe/Color:	White/Blanc/Bianco/Weiß/Blanco	
Density/Densité/Densita/Dichte/Densidad:	10.6 g/cm <sup>3</sup>	112 dwt/in <sup>3</sup>
Hardness/Dureté/Durezza/Härte/Dureza:	425	HV
Proof stress/Limite Élastique/Limite Elastico/ Dehngrenze/Limite de Elasticidad (0.2%):	900 MPa	130,500 psi
Elongation/Allongement/Allungamento/ Bruchdehnung/Alargamiento:	23	%
Melting range/Intervalle de fusion/Intervallo di fusione/ Schmelzintervall/Intervalo de fusion:	1100–1190 °C 2012–2174 °F	
Elastic Modulus/Module Élastique/Modulo Elastico/ Elastischer Modul/Módulo Elástico:	100 GPa	14,400 kpsi
Other components/Autres composants/Altri componenti/Weitere Bestandteile/ Otros componentes:	10.0 % Cu, 9.0 % Ga, 0.1 % Ir, 0.05 % B	

## Option® – Directions For Use

GB

### Indications: Dental Ceramic Alloy, Type 4

**Contraindications:** Do not use in patients with known hypersensitivity to one or several metals contained in the alloy.

**Warnings/Precautions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Use proper safety equipment and a certified industrially ventilated work area when processing this alloy.

**Adverse Reactions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

### Instructions:

**Investment:** Phosphate (carbon free).

**Burn out:** 700 °C (1300 °F), 30–120 min.

**Cast at:** 1290 °C (2354 °F), do not use carbon flux or crucibles.

**Finish with:** Carbide burs or ceramic bound stones.

**Pre-Clean:** Blast with alumina, 50–125 µm/< 3 bar (44 psi), steam clean or ultrasonic with distilled water.

**Oxidation:** 650 °C (1200 °F) to 1010 °C (1850 °F); no vacuum, hold 5 minutes.

**Post-Clean:** Do not blast, steam clean or ultrasonic with distilled water.

**Heat treatment:** None.

**Recommended Porcelain:** Ceramco® 3, Ceramco® II, Finesse®, Duceram® Kiss, Duceram® Plus.

**Pre-Solder:** Option® Pre/YPG; Flux: Fluoride.

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