MSDS

**R<sub>x</sub> only****DENTSPLY**  
CERAMCO

# Bioney<sup>®</sup>

Dental Ceramic Alloy, Type 4/Alliage céramique dentaire, Type 4/  
Leghe ceramiche dentale, Tipo 4/Dental-Keramiklegierung, Typ 4/  
Aleación cerámica dental, Tipo 4

90.0 % Au, 7.3 % Pt, 1.0 % Ag

CTE/CET/WAK 14.8 (500 °C), 15.1 (600 °C)

**REF 1810 0035****1 T.oz/31.1 g**

DENTSPLY CERAMCO  
DENTSPLY International Inc.  
570 West College Ave.  
York, PA 17405-0872  
1-800-487-0100  
1-717-699-4190

Color/Couleur/Colore/Farbe/Color:	Yellow/Jaune/Giallo/Gelb/Amarillo	
Density/Densité/Densita/Dichte/Densidad:	18.9 g/cm <sup>3</sup>	199 dwt/in <sup>3</sup>
Hardness/Dureté/Durezza/Härte/Dureza:	180	HV
Proof stress/Limite Élastique/Limite Elastico/ Dehngrenze/Limite de Elasticidad (0.2 %):	480 MPa	69,600 psi
Elongation/Allongement/Allungamento/ Bruchdehnung/Alargamiento:	9	%
Melting range/Intervalle de fusion/Intervallo di fusione/ Schmelzintervall/Intervalo de fusion:	1030–1095 °C 1890–2000 °F	
Elastic Modulus/Module Élastique/Modulo Elastico/ Elastischer Modul/Módulo Elástico:	76 GPa	11,000 kpsi
Other components/Autres composants/Altri componenti/Weitere Bestandteile/ Otros componentes:	0.5 % Zn, 0.5 % In, 0.7 % Fe, 0.05 % Re	

## Bioney® – Directions For Use

GB

### Indications: Dental Ceramic Alloy, Type 4

**Contraindications:** Do not use in patients with known hypersensitivity to one or several metals contained in the alloy.

**Warnings/Precautions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Use proper safety equipment and a certified industrially ventilated work area when processing this alloy.

**Adverse Reactions:** Exposure to alloy dust or fumes may cause eye irritation and/or respiratory complications. Reactions that may occur are allergic reactions against metals contained in the alloy or paraesthesia caused by electrochemical reactions. Systemic side effects caused by metals contained in the alloy have been reported in isolated cases.

### Instructions:

**Investment:** Phosphate.

**Burn out:** 700 °C (1300 °F), 30–120 min.

**Cast at:** 1190 °C (2175 °F).

**Finish with:** Carbide burs.

**Pre-Clean:** Blast with alumina, 50–125 µm/< 3 bar (44 psi), steam clean or ultrasonic with distilled water.

**Oxidation:** 650 °C (1200 °F) to 955 °C (1750 °F); no vacuum, no hold.

**Post-Clean:** Do not blast, steam clean or ultrasonic with distilled water.

**Recommended Porcelain:** Ceramco® 3, Ceramco® II, Ceramco® II Silver, Duceram® Kiss, Duceram® Plus.

**Heat treatment:** Harden following porcelain firing at 500 °C (930 °F) for 15 minutes, slow cool.

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