# silaplast

**FUTUR**

**GB**

Instructions for use

**Precision impression material, silicone based, condensation curing, very high consistency – putty**

Fig. 1

Fig. 2


Fig. 3

Fig. 4

## 1. Tray preparation

Select a tray type suited for desired application. For optimal adhesion of impression material, we recommend use of **sili adhesive**. Apply a thin layer of adhesive film to tray, then let dry for 60 seconds (Fig. 1).

## 2. Mixing

Fill measuring scoop evenly to top with **silaplast FUTUR** (Fig. 2). Smooth **silaplast FUTUR** onto mixing pad with a spatula and make a series of crosswise grooves (Fig. 3). Add recommended amount of catalyst (Fig. 4). Mix thoroughly, first using spatula and then by kneading between fingers until uniform colour is achieved (approx. 30 sec.).

## 3. Cleaning and disinfecting

After removing tray from mouth, rinse impression with lukewarm water. If desired, impression can then be disinfected by immersing into 2% glutaraldehyde for 15 min.

## 4. Pouring the model

We recommend waiting at least 30 minutes, but no more than 72 hours before pouring the model. All class III and IV dental plasters and standard model-casting plastics can be used.

## 5. Electro plating

The impression can be electro plated using all commercially available copper or silver plating baths.

## 6. Cleaning the tray

After setting, material can be dislodged with a blunt instrument. Remaining adhesive film can easily be removed with **sili Solvent** spray (FCHC-free; use only in well-ventilated rooms). Tray can then be cleaned and disinfected in the usual manner.

## Important hints


- We recommend using **silasoft® Normal**, **silasoft® Special** or **silasoft® direct** for corrective impression after making a preliminary impression with **silaplast FUTUR**.
- Overdosing of catalyst and increased temperatures accelerate, underdosing and decreased temperatures retard the setting time.
- Catalyst bottles or tubes should be closed tightly immediately after use.
- Cured impression materials are chemically inert – spots on clothing should be avoided.
- Avoid contact of catalyst with eyes, may cause irritation. In case of accidental contact rinse eyes with plenty of water and seek medical care immediately.

Further information:


Silicone based materials are proven a million times. On condition of a proper application, **undesired effects** are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For dental use only.

## silaplast FUTUR / cat<sup>p</sup> universal FUTUR

## silaplast FUTUR / cat<sup>f</sup> silaplast FUTUR

## Indications for use:

- Preliminary impressions for two-stage putty and wash techniques
- Situation impressions
- Check bites

## Technical data:

DIN EN ISO 4823 - Type 0

- **Mixing ratio:**  
**silaplast FUTUR / cat<sup>f</sup> silaplast FUTUR**  
1 scoop = 13.5 g /  
1 scoop = 0.50 g  
**silaplast FUTUR / cat<sup>p</sup> universal FUTUR**  
1 scoop = 13.5 g /  
6 cm paste = 0.30 g
- **Colour code:**  
**silaplast FUTUR:** white  
**cat<sup>p</sup> universal FUTUR:** blue  
**cat<sup>f</sup> silaplast FUTUR:** light blue
- **Mixing time:**  
approx. 30 sec.
- **Working time:**  
approx. 1 min. 15 sec.\*
- **Setting time:**  
**cat<sup>p</sup> universal FUTUR:**  
approx. 4 min. 45 sec.\*  
**cat<sup>f</sup> silaplast FUTUR:**  
approx. 3 min. 45 sec.\*
- **Time in the mouth:**  
**cat<sup>p</sup> universal FUTUR:**  
approx. 3 min. 30 sec.  
**cat<sup>f</sup> silaplast FUTUR:**  
approx. 2 min. 30 sec.
- **Strain in compression:**  
3.8 %
- **Recovery from deformation:** 98.60 %
- **Linear dimensional change:** 0.30 %
- **Application:**  
At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**  
At 15-25 °C / 59-77 °F,  
< 50 % rel. humidity

\* from initiation of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard these times.

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