silasoft® Special

GB
Instructions for use

Precision impression material, silicone based, condensation curing, low consistency – light bodied

Fig. 1


Fig. 2

Fig. 3

Fig. 4

1. Preparation of the impression

Depending on the impression technique select an appropriate tray. For optimal adhesion of the impression to the tray apply a thin layer of **Sili Adhesive Spray** or **Sili Adhesive** liquid onto the tray and let dry for 60 sec (fig. 1). Take preliminary impression with **silaplast FUTUR**. To facilitate the flowing out of excessive light bodied material from the impression cutting of escape vents and removal of interdental septi in the impression with the use of the special instrument **Deta-Cut** is recommended.

2. Mixing

Extrude **Silasoft® S** and **catp universal FUTUR** in equal lengths and thickness onto mixing pad (fig. 2). With the **DOSIMAT** system mixing of **silasoft® S** and catalyst paste is always automatically correct. Mix according to dosing ratio for approx. 30 sec. with spatula until homogeneous uniform colouring is achieved. Load preliminary impression. A normal syringe or disposal syringe can be used for bubble free application (fig. 3 + 4).

3. Disinfection

After removal rinse impression under lukewarm water. Disinfection for 15 min. in 2 % glutaraldehyde.

4. Model casting

Casting between 30 minutes and 24 hours after removal of the impression. Recommended materials: Dental plasters of class III and IV-also commercially available model acrylics.

5. Galvanisation

Impressions can be electroplated with commercially available copper and silver plating baths.

6. Cleaning of trays

Remove set material with a blunt instrument. **Sili Spray** or **Sili liquid** film can be easily removed with **Sili Spray Solvent** (free of FCHC). Work in well ventilated rooms. Clean und disinfect as usual.

Important working hints


- Increased temperatures accelerate, decreased temperatures retard the setting time.
- Catalyst tubes should be closed tightly immediately after use.
- Don't leave any residual material in the patient's mouth.
- Avoid staining of clothing. Hardened impression materials cannot be removed.
- Avoid contact of catalyst with eyes, may cause irritation. If necessary flush eyes with plenty of water and seek immediate medical care.

Other information:

Silicone impression materials are proven a million times over, unwanted effects cannot be expected in proper use. However, immune reactions, e.g. allergies or irritations cannot be excluded basically. In the case of doubt, we recommend performing an allergy test before application.

For use by trained specialists.

silasoft® S / catp universal FUTUR:

Indications for use:

- Corrective ("Wash")-material for double impression
- Reline impression

Technical Data:

DIN EN ISO 4823 – Type 3

- **Mixing ratio:**
silasoft® S / catp universal FUTUR
10 cm silasoft® S = 4.5 g / 10 cm catp = 0.5 g
- **Colour code:**
silasoft® S: rose-transparent
catp universal FUTUR: blue
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 1 min. 30 sec*
- **Setting time:** approx. 4 min. *
- **Time in mouth:** approx. 2 min. 30 sec.
- **Strain in compression:** 11.2 %
- **Recovery from deformation:** 98.5 %
- **Linear dimensional change:** 0.5 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity
- **Storage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

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