blue eco scan



GB
Instructions for use

Scannable laboratory putty, addition curing, smooth, non-sticky consistency

For the powderless 3D data acquisition of the antagonists for the CAD/CAM/CIM technologies (e.g. strip light projection, laser triangulation and CEREC procedures), excellent resolution for artifact-free image reproduction with high brightness and strong contrast. Easy mixing due to contrast colours and smooth consistency. Can be cut and milled very well, no breaking out of thin contact point areas.

Mixing and dispensing

Remove components A (white) and B (grey) with the corresponding measuring scoop and mix quickly in the ratio 1:1 (approx. 30 seconds), until a homogeneous mixture or uniform colour is achieved.

Important working hints

- For further work on the registration, observe the instructions of the scan system manufacturer.
- Clean and dry the registration before scanning.
- Only for dental **extraoral** use by qualified personnel.
- Do not bring the components A + B into contact before mixing, do not interchange measuring scoops and jar lids.
- Do not remove material by hand, do not wear latex gloves.

1 min. 30 sec.

Working time

4 min. 30 sec.

Setting time

Indications for use:

- Occlusal registration for powderless 3D data acquisition
- Conventional indications for laboratory putty
- Fabrication of small model segments

Technical Data:

- **Mixed volume:** 2 x 700 g (jars)
- **Mixing ratio:** 1:1
- **Colour code:** base: grey catalyst: white
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 1 min. 30 sec.*
- **Setting time:** approx. 4 min. 30 sec.*
- **Hardness Shore A:** 85
- **Linear dimensional change:** < 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering information:

blue eco scan
Standard packing **03114**
base, jar of 700 g
catalyst, jar of 700 g
2 measuring scoops

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