


ODUS CROWN FORMS MOULD CHART

In order to avoid confusion between ODUS-PELLA, ODUS-ACRYL, right and left crown forms, they are supplied filled with coloured cotton-wool:

white = PELLA right crowns
pink = PELLA left crowns
blue = ACRYL right crowns
yellow = ACRYL left crowns

Assortments supplied in:
Plexi-Boxes of 24 crowns (12 pairs)
Plexi-Boxes of 78 crowns (39 pairs)

Refill boxes of anyone number:
Boxes of 12 crowns (6 pairs)
Boxes of 24 crowns (12 pairs)
Boxes of 36 crowns (18 pairs)
Boxes of 60 crowns (30 pairs)



Numbering-System:

The numbers of all

Centrals _____ end with a 1
 Laterals _____ end with a 2
 Canines _____ end with a 3
 1st bicuspids _____ end with a 4
 2nd bicuspids _____ end with a 5
 upper molars _____ end with a 6
 lower molars _____ end with a 7